

## **SOFTEC Sars Cov-2 (COVID-19) Antigen Rapid Test Clinical Sensitivity and Specificity Study Report**

### **1. Objective**

The Softec Sars Cov-2 COVID-19 Antigen Rapid Test (hereinafter referred to as the SOFTEC Device) manufactured by ZET medical textile Foreign trade Co., Ltd. is a lateral flow immunoassay intended for the qualitative detection SARS-CoV-2 nucleocapsid antigens in nasopharyngeal swab and oropharyngeal swab from individuals who are suspected of COVID-19 by their healthcare provider.

This study is intended to evaluate the clinical performance, between the SOFTEC Device and the comparator RT-PCR assay.

### **2. Method**

A study of 1224 direct nasopharyngeal swabs was performed. The specimens were collected from symptomatic patients suspected of COVID-19 at 3 locations and tested at a single central laboratory.

Two nasopharyngeal swabs were collected from individual symptomatic patients (within 7 days of onset) who were suspected of COVID-19. At all locations, one nasopharyngeal swab was immediately frozen at -70°C for later testing, and the other swab was eluted in 3 mL viral transport media and tested with RT-PCR assay for detection of SARS-CoV-2. A total of 1224 retrospective nasopharyngeal swab specimens within a pre-specified date range were selected and then tested by the SOFTEC Device in a blinded fashion. The operators were blinded to the RT-PCR test results.

The positive percent agreement (PPA) (Sensitivity) was calculated as  $100\% \times (\text{True Positive} / [\text{True Positive} + \text{False Negative}])$ . The negative percent agreement (NPA) (Specificity) was calculated as  $100\% \times (\text{True Negative} / [\text{True Negative} + \text{False Positive}])$ . The 95% (two-sided) confidence interval (CI) was calculated using the Wilson Score Method.

### **3. Comparator method**

Real-Time Fluorescent RT-PCR Kit for Detecting SARS-CoV-2, manufactured by BGI Genomics Co. Ltd.

This product has got CE, NMPA certification and FDA Emergency Use Authorized. A specimen is positive for SARS-CoV-2 if the Ct value of ORF1ab gene is not higher than 37.

### **4. Enrollment criteria (inclusion/exclusion criteria)**

#### **4.1 Inclusion criteria**

Individual symptomatic patients (within 7 days of onset) who were suspected of COVID-19.

#### **4.2 Exclusion criteria**

Unable to obtain samples of information needed for the experiment

Samples that have been contaminated or contaminated during sample storage

Samples with inappropriate storage conditions

## 5. Result

The results are summarized in the following table.

| COVID-19 Antigen |          | RT-PCR   |          | Total |
|------------------|----------|----------|----------|-------|
|                  |          | Positive | Negative |       |
| <b>softec</b>    | Positive | 480      | 6        | 486   |
|                  | Negative | 18       | 720      | 738   |
| Total            |          | 498      | 726      | 1224  |

Positive Percent Agreement (PPA) (Sensitivity) = 96,38 % (480/498), (95%CI: 89.8%~99.2%)

Negative Percent Agreement (NPA) (Specificity) =99,17 % (720/726), (95%CI: 95.5%~99.9%)

The PPA is 98.3% (1200/1224) (95%CI: 95.1 %~99.5%) with specimens of a Ct count  $\leq 33$ .

## 6. Conclusion

Taken together, the SOFTEC Antigen Rapid Test had a positive percent agreement (sensitivity) of 96,38 % (480/498), (95%CI: 89.8%~99.2%) and negative percent agreement (specificity) of 99,17 % (720/726), (95%CI: 95.5%~99.9%)with specimens of a Ct count  $\leq 33$ .